

Standard	Dilution	Contact time	Temperature	Strains
EN 1276 - Bactericidal	ready to use	1 minute	20°C	<i>Pseudomonas Aeruginosa, Staphylococcus Aureus, Enterococcus Hirae, Escherichia Coli</i>
EN 1500 - Bactericidal	ready to use	1 minute	20°C	<i>Escherichia Coli</i>
EN 1650 - Yeasticidal	ready to use	1 minute	20°C	<i>Candida Albicans</i>
EN 13624 - Yeasticidal	ready to use	1 minute	20°C	<i>Candida Albicans</i>
EN 13727 - Bactericidal	ready to use	1 minute	20°C	<i>Pseudomonas Aeruginosa, Staphylococcus Aureus, Enterococcus Hirae, Escherichia Coli</i>
EN 14476 - Virucidal	ready to use	1 minute	20°C	<i>Murine Norovirus, Poliovirus Type 1, Adenovirus Type 5</i>
EN 12791 - Bactericidal	ready to use	3 minutes	20°C	<i>Resident skin microbial flora</i>

Material to be used exclusively by the Sutter Professional sales force